## H & A Canada Inc.

Salmonella

1160 Tapscott Road, Toronto, Ontario, Canada M1X 1E9 Tel: (416) 412-9518 Fax: (416) 293-9066 sales@hacanada.com



CP2005

## **Product Data**

Resveratrol 20%, 50%, 98%

| Product Information |                              |  |
|---------------------|------------------------------|--|
| Product Name        | Polygonum Cuspidatum Extract |  |
| Botanical Source    | Polygonum cuspidatum Sieb.   |  |
| Part Used           | Root                         |  |
| Solvent Used        | Ethanol & Water              |  |

| Analysis   | Specification                   | Test Method |
|------------|---------------------------------|-------------|
| Assay      | NLT 20%, 50% or 98% Resveratrol | HPLC        |
| Appearance | Brown yellow fine powder        | Visual      |

| Odor                   | Characteristic    | Organoleptic         |
|------------------------|-------------------|----------------------|
| Taste                  | Characteristic    | Organoleptic         |
| Sieve Analysis         | 100% pass 80 mesh | 80 Mesh Screen       |
| Bulk Density           | 45~70 g/100ml     | /                    |
| Loss on Drying         | ≤5.0%             | 5 g/105℃/2hrs        |
| Ash                    | ≤5.0%             | <b>2 g/525℃/3hrs</b> |
| Heavy Metals           | ≤10ppm            | Colorimetry          |
| Pb                     | <0.2ppm           | AFS                  |
| As                     | <1ppm             | AFS                  |
| Cd                     | <1ppm             | AFS                  |
| Hg                     | <1ppm             | AFS                  |
| Total Plate Count      | 1000 cfu/g Max    | CP2005               |
| Microbiological Contro | I                 |                      |
| Yeast & Mold           | 100 cfu/g Max     | CP2005               |
| E. Coli                | Negative          | CP2005               |

| <b>General Information</b> |   |
|----------------------------|---|
| Packing                    | Packed in 25 kg drums or 5 kg aluminium bags (98%)                |
| Storage                    | Store in a cool, dry place. Keep away from strong light and heat. |
| Sterilization              | This product is non IR and non ETO                                |

Negative

The information above is believed to be accurate and reliable. But no warranty or representation of any kind is made by us for its use. This information does not imply a recommendation to use in whatever way that may damage any products or infringe any existing patents or proprietary rights. In no way is H & A Canada Inc. liable or accountable for any claims or losses of any client or third party as a direct or indirect result of the use of this information.